

ELOQUENCE Hyaluronic Acid Dermal Fillers: Long-Lasting Efficacy

ZEEP Technology

Harnessing the power of Zero Endotoxin & BDDE Entire Process, ELOQUENCE ensures an endotoxin level below 0.1 EU/ml. By eliminating any detectable residual BDDE, we minimize potential side effects post-procedure.

The Downing Process

ELOQUENCE champions a unique Downing Process, producing a uniform particle size for smoother application. With this, professionals can achieve precision without particle agglomeration.

Safety First with Minimized BDDE Input

The Degree of Modification (MoD) is crucial for safety validation. Defined as 'MoD = full cross-linked BDDE + Pendant type BDDE', it gauges the modification ratio of HA by measuring the total BDDE in the HA Gel. ELOQUENCE 's MoD is 1/6 lower than other mono-type fillers, elevating its safety profile.

Maximized Cross-Linking Efficiency

Through an advanced physical cross-link process, ELOQUENCE attains optimal cross-linking. This is achieved while notably reducing the BDDE amount used as a chemical catalyst, resulting in a gel rich in viscoelastic properties.

Superior Viscoelasticity & Cohesiveness

When juxtaposed with leading filler products, ELOQUENCE stands tall. Its viscosity and cohesiveness are on par with U.S. products known for higher BDDE input, promising exceptional clinical outcomes.

Clinical Results & Efficacy Assessment

Independent expert panels and blinded investigators, using the WSRS results (score reductions compared to the baseline), have confirmed the noticeable



ELOQUENCE DEEP with HA and Lidocaine

For treating deep wrinkles, nasolabial folds and perfect for augmenting cheeks, forehead, and lips.



HA Concentration 24mg/ml



Lidocaine HCl 0.3%



Volume 1.1ml




2 sterile needles (27G x 2)



Duration of up to 12 months



CE certified



ELOQUENCE SUB-Q with HA and Lidocaine

For creating the most perfect augmentation of the cheekbones, nose, chin, and jawline.



HA Concentration 24mg/ml



Lidocaine HCl 0.3%



Volume 1.1ml




2 sterile needles (27G, 30G)



Duration of up to 12 months



CE certified



NEOM

A e s t h e t i c

ALL-IN-ONE LINE

Shape the New You,
Effortlessly and Easily



For Body and Face

ELOQUENCE SHAPE

ELOQUENCE CONTOUR



Easy to use, Easy to extract

ELOQUENCE SHAPE & ELOQUENCE CONTOUR

These two are all-in-one dermal fillers with high viscoelasticity that can be applied to both the face and body, which **offers the best value for both the practitioner and patient.**

Experience the ease and convenience of ELOQUENCE SHAPE: We've improved our formula to make it **easier to use, addressing any previous difficulties in extraction.**

With 95% of users expressing satisfaction with the quality of our HA fillers, 92% finding them easy to integrate into the dermis, and 91% reporting a smooth and uniform application, ELOQUENCE SHAPE & CONTOUR stand out as the go-to solution for effortless body contouring.

ELOQUENCE SHAPE

HA Concentration 24mg/ml

Shelf life : 3 years

Package : 1 Syringe 50ml

ELOQUENCE CONTOUR

HA Concentration 24mg/ml

Shelf life : 2 years

Package : 1 Syringe 10ml

LIDOCAINE 0.3%

EQ ELOQUENCE

