



NUCLEOFILL STRONG



PROMOITALIA
MEDICAL AESTHETICS

Features:

- 2.5% polynucleotide gel
- Molecular weight 2.500 KD
- 1.5 ml glass vial
- pH between 6.8 and 7.5
- Viscosity 48 - 65 Pa * s
- Two 30G 12mm needles included in the package

NUCLEOFILL STRONG



Indications

- ▶ Biorestructuring of the deep dermis, adipose tissue and bone tissue
- ▶ Prolonged moisturizing effect
- ▶ Intense lifting effect
- ▶ Treatment of small depressed scars and surgical scars

Nucleofill Strong is suitable for curative treatment in case of evident aging, when the skin has an accentuated roughness and a very loose tone.